

## SPECIFICATION

## Bis[(+)-tartrato]diantimonate(III) Dipotassium Trihydrate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (C <sub>8</sub> H <sub>4</sub> K <sub>2</sub> O <sub>12</sub> Sb <sub>2</sub> ·3H <sub>2</sub> O) Solubility in water pH of a 50g/l solution at 25°C Chloride (Cl) Sulfate (SO <sub>4</sub> ) Copper (Cu) Calcium (Ca) Zinc (Zn) Lead (Pb) Arsenic (As) Iron (Fe)	White, crystals ~ crystalline powder min.99.5%(mass/mass) to pass test 3.8~4.5 max.0.001%(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass) max.0.005%(mass/mass) max.5ppm(mass/mass) max.5ppm(mass/mass)

JIS K 8533:2012

(\*:Additional test performed by Wako but not required by JIS)