

Code No.02-1454

SPECIFICATION

1-Butanol

for Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20℃)	0.808~0.812g/mL
Refractive index 20℃	1.398~1.400
Absorbance (230nm)	max.0.20
Absorbance (240nm)	max.0.10
Absorbance (254nm)	max.0.03
Absorbance (270~400nm)	max.0.01
Water	max.0.1%
Fluorescence test	to pass test
Assay (capillary GC)	min.99.7%