

## SPECIFICATION

t-Butyl Methyl Ether, Super Dehydrated  
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	0.739~0.744g/mL
Water	max.0.001%
Residue after evaporation	max.10ppm
Acidity (as CH <sub>3</sub> COOH)	max.0.001%
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.5ppm
Assay (capillary GC)	min.99.5%