

Code No.03-0218

## SPECIFICATION

## Sodium p-Toluenesulfonchloramide Trihydrate

Guaranteed Reagent

| REQUIREMENT  | SPECIFICATION  |
|--|--|
| * Appearance   | White ~ slightly yellow, crystals ~ crystalline powder |
| Assay (C <sub>7</sub> H <sub>7</sub> ClNaO <sub>2</sub> S · 3H <sub>2</sub> O) | 98.0~104.0%(mass/mass)                                 |
| Solubility in water  | to pass test   |
| pH of a 50g/L solution at 25°C   | 8.0~11.0   |
| Substances darkened by sulfuric acid   | to pass test   |

JIS K 8318:2023(\*:Additional test performed by Wako but not required by JIS)