

Code No.03-0224

# SPECIFICATION

## Chloroacetic Acid

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals or mass
Solubility in water	to pass test
Melting point	61~64°C
Loss on drying over silica gel	under investigation
Chloride (Cl)	max.0.005%
Sulfate (SO <sub>4</sub> )	max.0.03%
Heavy metals (as Pb)	max.0.002%
Iron (Fe)	max.0.005%
Substances darkened by sulfuric acid	to pass test
Assay (after drying)	min.97.0%