

Code No.03-0399

SPECIFICATION

Copper, Powder
Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Reddish brown powder
Solubility in dil.nitric acid	to pass test
Sulfur compounds (as SO4)	max.0.02%
Silver (Ag)	max.0.005%
Lead (Pb)	max.0.01%
Arsenic (As)	max.0.002%
Antimony (Sb)	max.0.01%
Iron (Fe)	max.0.005%
Assay	min.99.5%