

Code No.03-0402

## SPECIFICATION

### Copper(II) Acetate Monohydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Bluish green, crystals ~ crystalline powder
Solubility in dil.acetic acid	to pass test
Chloride (Cl)	max.0.003%
Sulfate (SO <sub>4</sub> )	max.0.01%
Nitrogen compounds (as N)	max.0.05%
Sodium (Na)	max.0.2%
Potassium (K)	max.0.01%
Calcium (Ca)	max.0.01%
Zinc (Zn)	max.0.002%
Lead (Pb)	max.0.005%
Iron (Fe)	max.0.002%
Nickel (Ni)	max.0.002%
Assay	min.99.0%