

Code No.03-0413

SPECIFICATION

Copper (II) Chloride Dihydrate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Bluish green, crystals ~ crystalline powder
Solubility in ethanol	to pass test
Solubility in dil. hydrochloric acid	to pass test
Sulfate (SO ₄)	max.0.015%
Nitrogen compounds (as N)	max.0.01%
Arsenic (As)	max.4 ppm
Iron (Fe)	max.0.01%
Assay	min.97.0%