

Code No.03-0459

## SPECIFICATION

## Creatinine

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystals ~ powder
Solubility in water	to pass test
Chloride (Cl)	max.0.05%
Sulfate (SO <sub>4</sub> )	max.0.005%
Heavy metals (as Pb)	max.0.001%
Iron (Fe)	max.0.001%
Ammonium (NH <sub>4</sub> )	max.0.02%
Assay	min.99.0%