

Code No.03-0490

## SPECIFICATION

## Cupferron

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ pale yellowish brown, crystals ~
Solubility in water	crystalline powder
Specific absorbance[E1%1cm(281nm)]	to pass test
Loss on drying	min.515
Residue after ignition (as sulfate)	max.1.0%(mass/mass)
Sensitivity	max.0.05%(mass/mass)
	to pass test

JIS K 8289:2012(\*:Additional test performed by Wako but not required by JIS)