

Code No.03-1060

SPECIFICATION

Clofibrate

for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	1.139~1.146g/mL
Refractive index 20°C	1.500~1.505
Water	max.0.2%
Acidity (as C ₁₀ H ₁₁ ClO ₃)	max.0.2%
p-Chlorophenol (C ₆ H ₄ OH)	max.0.003%
Assay (capillary GC)	min.98.0%