

Code No.03-1872

## SPECIFICATION

Calcium Stearate, 2.0  $\mu$ m

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, powder
Particle diameter (d50)	max.2.0 $\mu$ m
Loss on drying at 105°C	max.5.0%
Acidity [as CH <sub>3</sub> (CH <sub>2</sub> ) <sub>16</sub> COOH]	max.0.5%
Assay (as Ca) (calculated on the dried basis)	6.0~7.0%
Assay (capillary GC)	under investigation