

## SPECIFICATION

## (E)-Chlorogenic Acid

for the Japanese Pharmacopoeia Crude Drugs Test (for Thin-Layer Chromatography)

REQUIREMENT	SPECIFICATION
Description (JP method)	White powder
Purity (JP method)	-
Related substances (JP method)	to pass test
* Solubility in methanol	to pass test
* Specific rotation [ $\alpha$ ] <sub>D20</sub> (c=0.4,H <sub>2</sub> O)	under investigation
* Water	under investigation
* Assay (HPLC)	min.98.0%

\*:Additional test performed by Wako