

Code No.03-2193

SPECIFICATION

Chloroform, Super Dehydrated, Amylene added for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Water	max.0.001%
Residue after evaporation	max.0.001%
Acidity (as HCl)	max.0.001%
Chloride (Cl)	max.1ppm
Free chlorine	to pass test (Cl:max.about 1ppm)
Substances darkened by sulfuric acid	to pass test
Assay (capillary GC)	min.99.5%
Amylene (stabilizer)	100~250ppm