

Code No.03-2395

SPECIFICATION

CultureSureR L-Tyrosine Disodium Salt Dihydrate, Animal-derived-free
for Cell Culture

REQUIREMENT	SPECIFICATION
Description	Light yellow, crystals~powder
Identification (1)	to pass test
Identification (2)	to pass test
Purity	-
Clarity of solution	to pass test
Water	12.7~14.7%
Assay	min.98.0%
Bacterial endotoxins	less than 10.0EU/g
Mycoplasma test	to pass test
Viable count test	max.100CFU/g