

Code No.03-2434

SPECIFICATION

CultureSureR Pyridoxine Hydrochloride, Animal-derived-free
for Cell Culture

REQUIREMENT	SPECIFICATION
Description	White to pale yellow, crystalline powder
Identification (1)	to pass test
Identification (2)	to pass test
Identification (3)	to pass test
pH	2.5~3.5
Purity	-
(1) Clarity and color of solution	to pass test
(2) Heavy metals	max.30ppm
(3) Related substances	to pass test
Loss on drying	max.0.30%
Residue on ignition	max.0.1%
Assay (after drying)	98.0~101.0%
Bacterial endotoxins	less than 50.0EU/g
Mycoplasma test	to pass test
Viable count test	max.50CFU/g