

Code No.03-2540

## SPECIFICATION

Cap B Solution [Acetic Anhydride-2,6-Lutidine-Acetonitrile (2:3:5)]  
for Nucleic Acid Synthesis

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Density (20°C)	0.883~0.889g/mL
Refractive index 20°C	1.398~1.402
Impurities	to pass test
Acetonitrile	49.0~51.0vol%
Acetic anhydride	19.0~21.0vol%
2,6-Lutidine	29.0~31.0vol%