

Code No.04-1935

## SPECIFICATION

## 2-(2-Butoxyethoxy)ethyl Acetate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ nearly colorless, clear liquid
Solubility in water	to pass test
Density (20°C)	0.975~0.981g/ml
Refractive index 20°C	1.424~1.428
Water	max.0.1%
Acidity (as CH <sub>3</sub> COOH)	max.0.01%
Assay (capillary GC)	min.98.0%