

Code No.04-2453

SPECIFICATION

Diethyl Ether

Infinity Pure

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.712~0.715g/ml
Water	max.0.1%
Residue after evaporation	max.2ppm
Acidity (as CH ₃ COOH)	max.10ppm
Aldehydes and ketones (as CH ₃ COCH ₃)	max.10ppm
Ethanol (C ₂ H ₅ OH)	max.0.05%
Substances darkened by sulfuric acid	to pass test
Assay (capillary GC)	min.99.5%