

Code No.04-2846

## SPECIFICATION

## Diethyl Ether for Pesticide Residue And Polychlorinated Biphenyl Analysis

Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
* Suitability for pesticide residue analysis (GC-FPD)	to pass test
Assay (C <sub>2</sub> H <sub>5</sub> OC <sub>2</sub> H <sub>5</sub> ) (GC)	min.99.5%(area/area)
Density (20°C)	0.712~0.714g/mL
Water	max.0.1%(mass/mass)
Residue after evaporation	max.2ppm(mass/mass)
Suitability for pesticide residue-PCB analysis	to pass test

JIS K 8357:2023(\*:Additional test performed by Wako but not required by JIS)