

SPECIFICATION

Dibasic Potassium Phosphate

Japanese Pharmaceutical Codex

REQUIREMENT	SPECIFICATION
* Description Identification pH Purity (1) Solubility (2) Chloride (3) Sulfate (4) Carbonate (5) Heavy metals (6) Sodium (7) Arsenic (8) Condensed phosphate Loss on drying Assay (after drying)	White, crystals or mass to pass test 8.6~9.3 - to pass test max.0.009% max.0.017% to pass test max.10ppm to pass test max.2ppm to pass test max.2.0% min.99.0%

Japanese Pharmaceutical Codex 2002(*:Reference item)