

Code No.04-2944

SPECIFICATION

Dibasic Sodium Phosphate Hydrate

Japanese Pharmacopoeia

REQUIREMENT	SPECIFICATION
* Description	Colorless or white crystals
Identification (1)	to pass test
Identification (2)	to pass test
Identification (3)	to pass test
pH	9.0~9.4
Purity	-
(1) Clarity and color of solution	to pass test
(2) Chloride	max.0.014%
(3) Sulfate	max.0.038%
(4) Carbonate	to pass test
(5) Heavy metals	max.10ppm
(6) Arsenic	max.2ppm
Loss on drying	57.0~61.0%
Assay (as anhydrous) (calculated on the dried basis)	min.98.0%

Japanese Pharmacopoeia(*:Reference item)