

Code No.04-3319

## SPECIFICATION

### Diisopropyl Ether, Super Dehydrated for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.720~0.730g/ml
Refractive index 20°C	1.366~1.371
Water	max.0.001%
Residue after evaporation	max.0.02%
Acidity (as CH <sub>3</sub> CH <sub>2</sub> COOH)	max.0.01%
Peroxide (as H <sub>2</sub> O <sub>2</sub> )	max.0.001%
Assay (capillary GC)	min.99.0%