

Code No.04-3457

# SPECIFICATION

## Dichloromethane-d<sub>2</sub>, 99.8% (High purity) for Nuclear Magnetic Resonance

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Water	max.0.005%
Degree of deuteration	min.99.8%
Assay (capillary GC)	min.99.9%
Impurities (NMR) (1H)	to pass test
Impurities (NMR) (13C)	to pass test
Impurities (NMR) (31P)	to pass test
Impurities (NMR) (19F)	to pass test