

Code No.05-0080

## SPECIFICATION

## Ethyl Crotonate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Solubility in ethanol	to pass test
Density (20°C)	0.917~0.924g/ml
Refractive index 20°C	1.422~1.428
Water	max.0.4%
Acidity (as CH <sub>3</sub> CH:CHCOOH)	max.0.2%
Assay (capillary GC)	min.95.0%