

# SPECIFICATION

## 2-Methoxyethanol for Amino Acid Automated Analysis

| REQUIREMENT                           | SPECIFICATION          |
|---------------------------------------|------------------------|
| Appearance                            | Colorless clear liquid |
| Solubility in sodium acetate solution | to pass test           |
| Density (20°C)                        | 0.962~0.967g/ml        |
| Water                                 | max.0.2%               |
| Acidity (as CH <sub>3</sub> COOH)     | max.0.005%             |
| Peroxides                             | to pass test           |
| Ninhydrin reactive substances         | to pass test           |
| Assay (capillary GC)                  | min.99.0%              |