

Code No.05-0110

## SPECIFICATION

## 2-Methoxyethanol

Wako 1st Grade

| REQUIREMENT                          | SPECIFICATION          |
|--------------------------------------|------------------------|
| Appearance                           | Colorless clear liquid |
| Solubility in water                  | to pass test           |
| Density (20°C)                       | 0.962~0.967g/mL        |
| Refractive index 20°C                | 1.401~1.404            |
| Water                                | max.0.1%               |
| Substances darkened by sulfuric acid | to pass test           |
| Assay (capillary GC)                 | min.98.0%              |