

Code No.05-0384

## SPECIFICATION

### Ethylene Glycol Diglycidyl Ether for Electron Microscopy

| REQUIREMENT         | SPECIFICATION                         |
|---------------------|---------------------------------------|
| Appearance          | Colorless ~ pale yellow, clear liquid |
| Viscosity (25°C)    | 15~25mPa·s                            |
| Solubility in water | to pass test                          |
| Solubility in NSA   | to pass test                          |
| Density (20°C)      | 1.14~1.16g/mL                         |
| Water               | max.0.3%                              |
| Practical test      | to pass test                          |
| Epoxy equivalent    | 110~120                               |