

SPECIFICATION

Ethanol(99.5) for Pesticide Residue And Polychlorinated Biphenyl Analysis
Concentration 300

| REQUIREMENT | SPECIFICATION |
|---|------------------------|
| * Appearance | Colorless clear liquid |
| * Suitability for pesticide residue analysis (GC-FPD) | to pass test |
| Assay (C ₂ H ₅ OH) (GC) | min.99.5%(area/area) |
| Density (20°C) | 0.789~0.791g/mL |
| Water | max.0.3%(mass/mass) |
| Residue after evaporation | max.2ppm(mass/mass) |
| Suitability for pesticide residue-PCB analysis | to pass test |

JIS K 8093:2023(*:Additional test performed by Wako but not required by JIS)