

Code No.05-0596

SPECIFICATION

Ethyl Propiolate

Wako Special Grade

| REQUIREMENT | SPECIFICATION |
|-----------------------|---------------------------------------|
| Appearance | Colorless ~ pale yellow, clear liquid |
| Density (20°C) | 0.963~0.968g/mL |
| Refractive index 20°C | 1.409~1.414 |
| Water | max.0.2% |
| Acidity (as HCCCCOOH) | max.0.3% |
| Assay (capillary GC) | min.97.0% |