

Code No.05-0632

SPECIFICATION

N-Ethyldiisopropylamine for Peptide Synthesis

| REQUIREMENT | SPECIFICATION |
|--------------------------------------------------------------------------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Water | max.0.01% |
| Ethylamine (CH ₃ CH ₂ NH ₂) | max.0.0015% |
| N-Ethylisopropylamine [(CH ₃) ₂ CHNHC ₂ H ₅] | max.0.015% |
| Isopropylamine [(CH ₃) ₂ CHNH ₂] | max.0.005% |
| Assay (capillary GC) | min.99.0% |