

SPECIFICATION

Formic Acid
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (HCOOH)	min.98.0%(mass/mass)
Solubility in water	to pass test
Residue after evaporation	max.0.002%(mass/mass)
Acetic acid (CH ₃ COOH) (GC)	max.0.3%(mass/mass)
Chloride (Cl)	max.5ppm(mass/mass)
Sulfate (SO ₄)	max.0.005%(mass/mass)
Sulfite	to pass test
Copper (Cu)	max.3ppm(mass/mass)
Lead (Pb)	max.3ppm(mass/mass)
Iron (Fe)	max.3ppm(mass/mass)

JIS K 8264:2020(*:Additional test performed by Wako but not required by JIS)The assay of this product is approximately 98.0%.