

## SPECIFICATION

## 5-Fluorouracil

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ powder
Absorption maximum	263~269nm
Loss on drying at 105°C	max.0.3%
Heavy metals (as Pb)	max.0.002%
Arsenic (As)	max.2ppm
5-Chlorouracil (C <sub>4</sub> H <sub>3</sub> ClN <sub>2</sub> O <sub>2</sub> )	max.0.1%
Assay (HPLC)	min.98.5%