

Code No.06-0472

SPECIFICATION

0.1vol% Formic Acid-Acetonitrile for LC/MS

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (200nm)	max.2.00
Absorbance (210nm)	max.2.50
Absorbance (220nm)	max.2.50
Absorbance (230nm)	max.1.50
Absorbance (240nm)	max.0.80
Absorbance (254nm)	max.0.05
Absorbance (260nm)	max.0.02
Absorbance (270~400nm)	max.0.01
Water	max.0.05%
Fluorescence test	to pass test
Suitability for LC/MS analysis	to pass test
Formic acid	0.14~0.18%