

SPECIFICATION

Florisol, 150~250 μ m (60~100mesh)
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Particle size (150~250 μ m)	min.60%
Particle size (coarser than 250 μ m)	max.40%
Particle size (finer than 150 μ m)	max.20%
Bulk density	under investigation
Loss on drying at 180°C	max.2.0%
pH (50g/l aqueous slurry,25°C)	8.6~10.2
Methanol dipping test	to pass test