

# SPECIFICATION

Florisil, 75~150  $\mu\text{m}$  (100~200mesh)  
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Particle size (75~150 $\mu\text{m}$ )	min.60%
Particle size (coarser than 150 $\mu\text{m}$ )	max.20%
Particle size (finer than 75 $\mu\text{m}$ )	max.40%
Bulk density	0.4~0.7g/mL
Loss on drying at 180°C	max.3.0%
pH(50g/L aqueous slurry, 25°C)	8.6~10.2
Methanol dipping test	to pass test