

Code No.06-0653

SPECIFICATION

Fumonisin Mixture Standard Solution [B1, B2 each 50 $\mu\text{g}/\text{mL}$ Acetonitrile:Water(1:1)
Solution]
for Mycotoxin Assay

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Practical test	to pass test
Fumonisin B1	49.2~50.8 $\mu\text{g}/\text{mL}$
Fumonisin B2	49.2~50.8 $\mu\text{g}/\text{mL}$