

SPECIFICATION

Glutaric Acid

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystalline powder
Solubility in water	to pass test
Melting point	95~99°C
Chloride (Cl)	max.0.001%
Sulfate (SO ₄)	max.0.003%
Copper (Cu)	max.0.001%
Lead (Pb)	max.0.001%
Iron (Fe)	max.0.001%
Assay (capillary GC)	under investigation
Assay (volumetric analysis)	min.98.0%