

Code No.07-0086

SPECIFICATION

Glycylglycine

Wako Special Grade

| REQUIREMENT | SPECIFICATION |
|-------------------------------------|--------------------------|
| Appearance | White crystalline powder |
| Solubility in water | to pass test |
| Loss on drying at 105°C | max.0.2% |
| Residue after ignition (as sulfate) | max.0.05% |
| Chloride (Cl) | max.0.005% |
| Sulfate (SO ₄) | max.0.005% |
| Heavy metals (as Pb) | max.0.001% |
| Iron (Fe) | max.0.003% |
| Thin layer chromatography test | to pass test |
| Assay | min.98.0% |