

SPECIFICATION

3-(o-Methoxyphenoxy)-1,2-propanediol

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystalline powder ~ powder
Solubility in water	to pass test
Melting point	80~85°C
Loss on drying at 60°C in vacuum	max.0.5%
Residue after ignition (as sulfate)	max.0.1%
pH of a 10g/l solution at 25°C	5.0~7.0
Guaiacol	to pass test
Analogous substances	to pass test
Assay	98.0~102.0%