

SPECIFICATION

[6]-Gingerol

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

REQUIREMENT	SPECIFICATION
Description (JP method)	Yellowish white ~ yellow liquid or solid
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
Unity of peak (JP method)	to pass test
Assay (qNMR) (JP method)	under investigation
* Solubility in methanol	to pass test
* Water	max.5.0%
* Assay (HPLC)	min.98.0%

This product is suited for [6]-Gingerol for assay (for assay 2) and thin-layer chromatography of general tests on The Japanese Pharmacopoeia. (*:Additional test performed by Wako)