

Code No.08-0015

SPECIFICATION

Heptane

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [CH ₃ (CH ₂) ₅ CH ₃] (GC) Density (20°C) Refractive index 20°C Residue after evaporation Substances darkened by sulfuric acid	Colorless clear liquid min.99.0%(mass/mass) 0.682~0.686g/ml 1.386~1.389 max.0.001%(mass/mass) to pass test

JIS K 9701:2013(*:Additional test performed by Wako but not required by JIS)