

Code No.08-0015

SPECIFICATION

Heptane
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay [$\text{CH}_3(\text{CH}_2)_5\text{CH}_3$] (GC)	min.99.0%(mass/mass)
Density (20°C)	0.682~0.686g/ml
Refractive index 20°C	1.386~1.389
Residue after evaporation	max.0.001%(mass/mass)
Substances darkened by sulfuric acid	to pass test

JIS K 9701:2013(*:Additional test performed by Wako but not required by JIS)