

Code No.08-0121

SPECIFICATION

Hydroquinone

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ nearly white, crystals
Assay	min.99.0%
Solubility in dil.acetic acid	to pass test
Melting point	171~174℃
Residue after ignition (as sulfate)	max.0.02%
Chlorine compounds (as Cl)	max.0.001%
Sulfate (SO ₄)	max.0.01%
Lead (Pb)	max.0.001%
Iron (Fe)	max.0.001%
Catechol (C ₆ H ₆ O ₂)	max.0.3%

JIS K 8738-1994

*:Additional test performed by Wako but not required by JIS