

## SPECIFICATION

## 8-Quinolinol

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay Solubility in ethanol Solubility in acetic acid Melting point Residue after ignition (as sulfate) Chloride (Cl) Sulfate (SO <sub>4</sub> ) Suitability for magnesium determination	White ~ pale yellow, crystals ~ crystalline powder min.99.0%  to pass test  to pass test  72~76°C  max.0.05%  max.0.001%  max.0.01%  to pass test

JIS K 8775-1994

\*:Additional test performed by Wako but not required by JIS