

Code No.08-0176

# SPECIFICATION

## Phosphinic Acid

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Phosphate (PO <sub>4</sub> )	max.0.1%
Sulfate (SO <sub>4</sub> )	max.0.1%
Heavy metals (as Pb)	max.0.002%
Barium (Ba)	max.0.02%
Arsenic (As)	max.0.001%
Oxalic acid	to pass test
Assay	30.0~32.0%