

Code No.08-0335

SPECIFICATION

Heptane

for Amino Acid Sequence Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (210nm)	max.0.80
Absorbance (230nm)	max.0.10
Absorbance (254~400nm)	max.0.01
Water	max.0.005%
Aldehydes (as CH ₃ CHO)	max.1ppm
Peroxides (as H ₂ O ₂)	max.1ppm
Assay (capillary GC)	min.99.0%