

Code No.08-0390

SPECIFICATION

o-Phenylphenol Standard

for Food Additives Test

| REQUIREMENT | SPECIFICATION |
|---|---|
| Appearance | White ~ slightly brown, crystals ~ powder |
| Solubility in ethanol | to pass test |
| Melting point | 55~59℃ |
| Water | max.0.2% |
| Biphenyl (C ₆ H ₅ C ₆ H ₅) | max.0.1% |
| Thin layer chromatography test | to pass test |
| Assay (capillary GC) | min.99.5% |