

Code No.08-0550

# SPECIFICATION

## Heptane

for SHIMADZU Protein Peptide Sequencer

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (210nm)	max.0.80
Absorbance (230nm)	max.0.10
Absorbance (254~400nm)	max.0.01
Water	max.0.005%
Aldehydes (as CH <sub>3</sub> CHO)	max.1ppm
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.1ppm
Suitability for amino acids sequence analysis	to pass test
Assay (capillary GC)	min.99.0%