

Code No.08-0744

SPECIFICATION

Hydrogen Peroxide

Super Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Residue after evaporation	max.10ppm
Chloride (Cl)	max.1ppm
Phosphate (PO ₄)	max.1ppm
Sulfate (SO ₄)	max.1ppm
Sodium (Na)	max.0.3ppm
Potassium (K)	max.0.3ppm
Copper (Cu)	max.0.003ppm
Calcium (Ca)	max.0.05ppm
Zinc (Zn)	max.0.01ppm
Cadmium (Cd)	max.0.001ppm
Mercury (Hg)	max.0.003ppm
Lead (Pb)	max.0.05ppm
Arsenic (As)	max.0.005ppm
Chromium (Cr)	max.0.01ppm
Manganese (Mn)	max.0.05ppm
Iron (Fe)	max.0.05ppm
Cobalt (Co)	max.0.003ppm
Nickel (Ni)	max.0.003ppm
Assay	30.0~35.5%